

A rapid test for the qualitative detection of Marijuana in human urine.
For professional in vitro diagnostic use only.

INTENDED USE

The THC Rapid Test Cassette (Urine) is a rapid chromatographic immunoassay for the detection of 11-nor- Δ^9 -THC-9 COOH (THC metabolite) in human urine at a cut-off concentration of 600ng/mL. This assay provides only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrophotometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are used.

SUMMARY

THC (Δ^9 -tetrahydrocannabinol) is the primary active ingredient in cannabinoids (Marijuana). When smoked or orally administered, it produces euphoric effects. Users have impaired short term memory and slowed learning. Users may also experience transient episodes of confusion and anxiety. Long term relatively heavy use may be associated with behavioral disorders. The peak effect of smoking Marijuana occurs in 20-30 minutes and the duration is 90-120 minutes after one cigarette. Elevated levels of urinary metabolites are found within hours of exposure and remain detectable for 3-10 days after smoking. The main metabolite excreted in the urine is 11-nor- Δ^9 -tetrahydrocannabinol-9-carboxylic acid (Δ^9 -THC-COOH). The THC Rapid Test Cassette (Urine) is a rapid urine screening test that can be performed without the use of an instrument. The test utilizes a monoclonal antibody to selectively detect elevated levels of Marijuana in urine. The THC Rapid Test Cassette (Urine) yields a positive result when the concentration of Marijuana in urine exceeds 600ng/mL. This is the suggested screening cut-off for positive specimens set by the Substance Abuse and Mental Health Services Administration (SAMHSA, USA).

PRINCIPLE

THC Rapid Test Cassette (Urine) is a rapid chromatographic immunoassay based on the principle of competitive binding. Drugs which may be present in the urine specimen compete against the drug conjugate for binding sites on the antibody. During testing, a urine specimen migrates upward by capillary action. Marijuana, if present in the urine specimen below 600ng/mL, will not saturate the binding sites of the antibody coated particles in the strip. The antibody coated particles will then be captured by immobilized THC conjugate and a visible colored line will show up in the test line region. The colored line will not form in the test line region if the Marijuana level is above 600ng/mL because it will saturate all the binding sites of anti-Marijuana antibodies. A drug-positive urine specimen will not generate a colored line in the test line region because of drug competition, while a drug-negative urine specimen or a specimen containing a drug concentration less than the cut-off will generate a line in the test line region. To serve as a procedural control, a colored line will always appear at the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

REAGENTS

The test contains mouse monoclonal anti-THC antibody-coupled particles and THC-protein conjugate. A rabbit antibody is employed in the control line system.

PRECAUTIONS

- For medical and other professional in vitro diagnostic use only. Do not use after the expiration date.
- The test should remain in the sealed pouch until use.
- All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The used test should be discarded according to local regulations.

STORAGE AND STABILITY

Store as packaged at room temperature or refrigerated (2-30°C). The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. **DO NOT FREEZE.** Do not use beyond the expiration date.

SPECIMEN COLLECTION AND PREPARATION

Urine Assay
The urine specimen must be collected in a clean and dry container. Urine collected at any time of the day may be used. Urine specimens exhibiting visible particles should be centrifuged, filtered, or allowed to settle to obtain clear specimen for testing.


Specimen Storage
Urine specimens may be stored at 2-8°C for up to 48 hours prior to testing. For long-term storage, specimens may be frozen and stored below -20°C. Frozen specimens should be thawed and mixed before testing.

MATERIALS

- Test Cassettes
- Package insert
- Droppers
- Specimen collection container
- Materials Required but Not Provided
- Timer

DIRECTIONS FOR USE

- Allow the test, urine specimen, and/or controls to reach room temperature (15-30°C) prior to testing.
- Bring the pouch to room temperature before opening it. Remove the test cassette from the sealed pouch and use it within one hour.
- Place the test cassette on a clean and level surface. Hold the dropper vertically and transfer 3 full drops of urine (approx. 120 μ L) to the specimen well (S) of the test cassette, and then start the timer. Avoid trapping air bubbles in the specimen well (S). See the illustration below.
- Wait for the color line(s) to appear. The result should be read at 5 minutes. Do not interpret the result after 10 minutes.



INTERPRETATION OF RESULTS

(Please refer to the illustration above)
NEGATIVE: Two lines appear. One color line should be in the control region (C), and another apparent color line should be in the test region (T). This negative result indicates that the Marijuana concentration is below the detectable level of 600ng/mL.
***NOTE:** The intensity of the color in the test line region (T) may vary depending on the concentration of 11-nor- Δ^9 -THC-9 COOH (THC metabolite) present in the specimen. Therefore, any shade of color in the test line region (T) should be considered negative.
POSITIVE: One color line appears in the control region (C). No line appears in the test region (T). This positive result indicates that the Marijuana concentration is above the detectable level of 600ng/mL.
INVALID: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new Test Cassette. If the problem persists, discontinue using the Test Cassette immediately and contact your local distributor.

QUALITY CONTROL

A procedural control is included in the test. A color line appearing in the control region (C) is considered an internal procedural control. It confirms sufficient specimen volume and correct procedural technique. Control standards are not supplied with this Test Cassette; however it is recommended that positive and negative controls be tested as good laboratory testing practices to confirm the test procedure and to verify proper test performance.

LIMITATIONS

- THC Rapid Test Cassette (Urine) provides only a qualitative, preliminary analytical result. A secondary analytical method must be used to obtain a confirmed result. Gas chromatography/mass spectrophotometry (GC/MS) is the preferred confirmatory method.^{1,2}
- It is possible that technical or procedural errors, as well as other interfering substances in the urine specimen may cause erroneous results.
- Adulterants, such as bleach and/or alum, in urine specimens may produce erroneous results regardless of the analytical method used. If adulteration is suspected, the test should be repeated with another urine specimen.
- A positive result indicates presence of the drug or its metabolites but does not indicate level of intoxication, administration route or concentration in urine.
- A negative result may not necessarily indicate drug-free urine. Negative results can be obtained when drug is present but below the cut-off level of the test.
- Test does not distinguish between drugs of abuse and certain medications.

EXPECTED VALUES

This negative result indicates that the Marijuana concentration is below the detectable level of 600ng/mL. Positive result means the concentration of Marijuana is above the level of 600ng/mL. The THC Rapid Test Cassette has a sensitivity of 600ng/mL.

PERFORMANCE CHARACTERISTICS

Accuracy
A side-by-side comparison was conducted using The THC Rapid Test Cassette and GC/MS at the cut-off of 600ng/mL. Testing was performed on 250 clinical specimens previously collected from subjects present for Drug Screen Testing. The following results were tabulated:

Method	Results	GC/MS		Total Results
		Positive	Negative	
The THC Rapid Test Cassette	Positive	81	5	86
	Negative	5	159	164
Total Results		86	164	250
% Agreement with this Rapid Test		94.2%	96.9%	96.0%

Analytical Sensitivity
A drug-free urine pool was spiked with 11-nor- Δ^9 -Tetrahydrocannabinol-9-COOH at the following concentrations: 0ng/mL, 300ng/mL, 450ng/mL, 600ng/mL, 750ng/mL, 900ng/mL and 1800ng/mL. The result demonstrates >99% accuracy at 50% above and 50% below the cut-off concentration. The data are summarized below:

11-nor- Δ^9 -THC-9 COOH Concentration (ng/mL)	Percent of Cut-off	n	Visual Result	
			Negative	Positive
0	0	30	30	0
300	-50%	30	30	0
450	-25%	30	27	4
600	Cut-off	30	14	16
750	+25%	30	4	26
900	+50%	30	0	30
1800	3X	30	0	30

The following table lists compounds and their respective concentrations in urine that yield a positive result in The THC Rapid Test Cassette (Urine) at 5 minutes.

Compound	Concentration (ng/mL)
Cannabinol	400,000
11-nor- Δ^8 - THC-9 COOH	400
11-nor- Δ^8 - THC-9 COOH	600
Δ^8 - THC	200,000
Δ^8 - THC	200,000

Precision

A study was conducted at three hospitals by laypersons using three different lots of product to demonstrate the within run, between run and between operator precision. An identical panel of coded specimens containing, according to GC/MS, no 11-nor- Δ^9 -Tetrahydrocannabinol-9-carboxylic acid, 25% 11-nor- Δ^9 -Tetrahydrocannabinol-9-carboxylic acid above and below the cut-off, and 50% 11-nor- Δ^9 -Tetrahydrocannabinol-9-carboxylic acid above and below the 600ng/mL cut-off was provided to each site. The following results were tabulated:

11-nor- Δ^9 -THC-9 COOH Concentration (ng/mL)	n per site	Site A		Site B		Site C	
		-	+	-	+	-	+
0	10	10	0	10	0	10	0
300	10	10	0	10	0	10	0
450	10	8	2	9	1	9	1
750	10	2	8	3	7	1	9
900	10	0	10	0	10	0	10

Effect of Urinary Specific Gravity
Fifteen urine specimens of normal, high, and low specific gravity ranges were spiked with 300ng/mL and 900ng/mL of 11-nor- Δ^9 -Tetrahydrocannabinol-9-carboxylic acid. The THC Rapid Test Cassette (Urine) was tested in duplicate using the fifteen neat and spiked urine specimens. The results demonstrate that varying ranges of urinary specific gravity do not affect the test results.

Effect of Urinary pH
The pH of an aliquoted negative urine pool was adjusted to a pH range of 5 to 9 in 1 pH unit increments and spiked with 11-nor- Δ^9 -Tetrahydrocannabinol-9-carboxylic acid to 300ng/mL and 900 ng/mL. The spiked, pH-adjusted urine was tested with the THC Rapid Test Cassette (Urine) in duplicate. The results demonstrate that varying ranges of pH do not interfere with the performance of the test.

Cross-Reactivity
A study was conducted to determine the cross-reactivity of the test with compounds in either drug-free urine or Marijuana positive urine. The following compounds show no cross-reactivity when tested with The THC Rapid Test Cassette (Urine) at a concentration of 100 μ g/ml.

Non Cross-Reacting Compounds		
4-Acetamidophenol	Deoxycorticosterone	(+) 3,4-Methylenedioxy-amphetamine
Acetophenetidin	Dextromethorphan	Prednisone
N-Acetylprocainamide	Diazepam	(+) 3,4-Methylenedioxy-methamphetamine
Acetylsalicylic acid	Diclofenac	Promazine
Aminopyrine	Difenhydramine	Promethazine
Amitypyline	Digoxin	Methyprylon
Amobarbital	Diphenhydramine	Morphine-3-
Amoxicillin	Doxylamine	D-glucuronide
Ampicillin	Egonine hydrochloride	D-Pseudoephedrine
L-Ascorbic acid	Egonine methylester	Naloxone
D,L-Amphetamine	(-)- ψ -Ephedrine	Naltrexone
L-Amphetamine	Erythromycin	Naproxen
Apomorphine	-Estradiol	Niacinamide
Aspartame	Estrone-3-sulfate	Nifedipine
Atropine	Ethyl-p-aminobenzoate	Norcodeine
Benzilic acid	Fenoprofen	Furosemide
Benzoic acid	Fenpropion	Gentic acid
Benzoylcegonine	Gentic acid	Hemoglobin
Benzphetamine	Hemoglobin	Hydrochlorothiazide
Bilirubin	Hydroalazine	Hydrocodone
(\pm)-Brompheniramine	Hydrochlorothiazide	Hydrocortisone
Caffeine	Hydrocodone	O-Hydroxyhippuric acid
Cannabidiol	Hydrocortisone	3-Hydroxytyramine
Chloralhydrate	Chlorpheniramine	Ibuprofen
Chloramphenicol	Chlorothiazide	Imipramine
Chloridiazepoxide	(\pm) Chlorpheniramine	Iproniazid
Chlorothiazide	Chlorpromazine	(\pm) - Isoproterenol
Cocaine hydrochloride	Chlorquine	Isoxsuprine
Cocaine	Cholesterol	Ketamine
Codine	Clomipramine	Ketoprofen
Cortisone	Clonidine	Labetalol
(-) Cotinine	Cocaine hydrochloride	Levorphanol
Creatinine	Codine	Loperamide
	Cortisone	Maprotiline
	(-) Cotinine	Meprobamate
	Creatinine	Methadone
		Methoxyphenamine
		Phenobarbital
		Phenylephrine
		-Phenylethylamine
		Phenylpropanolamine
		Prednisolon
		Procaine
		Promazine
		Promethazine
		D,L-Propranolol
		D-Proxophene
		D-Pseudoephedrine
		Quinidine
		Quinine
		Ranitidine
		Salicylic acid
		Secobarbital
		Serotonin (5-Hydroxytryptamine)
		Sulfamethazine
		Sulindac
		Temazepam
		Tetracycline
		Tetrahydrocortisone,
		3-Acetate
		Tetrahydrocortisone
		3-(D-glucuronide)
		Tetrahydrozoline
		Thebaine
		Thiamine
		Thioridazine
		D, L-Thyroxine
		Tolbutamide
		Triamterene
		Trifluoperazine
		Trimethoprim
		Trimipramine
		Triptamine
		D, L-Tryptophan
		Tyramine
		D, L-Tyrosine
		Uric acid
		Verapamil
		Zomepirac

BIBLIOGRAPHY

- Hawks RL, CN Chiang. Urine Testing for Drugs of Abuse. National Institute for Drug Abuse (NIDA), Research Monograph 73, 1986
- Baselt RC. Disposition of Toxic Drugs and Chemicals in Man. 2nd Ed. Biomedical Publ., Davis, CA. 1982; 488

Index of Symbols	
	Attention, see instructions for use
	For in vitro diagnostic use only
	Store between 2-30°C
	Do not use if package is damaged
	Tests per kit
	Use by
	Lot Number
	Authorized Representative
	Do not reuse
	Catalog #

ACRO BIOTECH, INC
9500 Seventh Street, Unit M, Rancho Cucamonga, CA 91730, U.S.A.
Tel.#: +1 (909) 466-6857 Fax #.: +1 (909) 466-6892
http://www.acrobiotech.com

CE **EC REP**
MedNet GmbH
Eckstrasse 10
48163 Muenster
Germany

DN: 304000
Rev.Date: 2016-06-21